



Calculations for Ion-Selective Measurements

1. Calculation of the Interference Error

The Interference Error is calculated from the:

The Ion Concentration of the Ion to be Measured,

The Ion Concentration of the Interfering Ion,

The Selectivity Coefficient of the Interference.

It will be best demonstrated by an Example:

Interference Error for the Interference of Potassium on an Ammonium Measurement:

Potassium Concentration in the sample: 50 ppm.

Selectivity Coefficient for Potassium on the Ammonium Electrode: 0.11

- a) Measuring 20 ppm Ammonium:
Error = $50 \times 0.11 / 20 = 0.275 = 27.5 \%$
- b) Measuring 100 ppm Ammonium:
Error = $50 \times 0.11 / 100 = 0,055 = 5.5 \%$

Fortunately many Selectivity Coefficients are much lower!

But some of them can be much higher, then the interfering Ion must be absent (Concentration less than 10×10^{-4} ppm).

2. Compensation of the Interference

2.1)

If the Concentration of the Interfering Ion in the sample is known and is constant, then one can do a Calibration with the Interfering Ion present in the Calibrating Solution.

2.2)

Other Methods are the so called "Incremental Methods":

- a) Standard Addition Method, where a known volume of a Calibration Standard is added to a known volume of Sample and the electrode signal before and after the addition is measured.
- b) Sample Addition Method, where a known volume of Sample is added to a known volume of a Calibration Standard and the electrode signal before and after the addition is measured.

2.3)

Making a Multi-Parameter Calibration and Measuring the Concentration of the Interfering Ion. Then using a Computer Algorithm to Calculate the Ion Concentration of the Ion to be measured, from both Electrode signals and the Multi-Parameter Calibration.